

BioSerenity | LR – PSG Data Compatibility Matrix

Updated 2022

Company	Product Name	Show list of tests for uploading	Find a sleep study on the local computer	Automatically extract patient & study info from PSG record	Download to the folder for native PSG	Extract scoring stats & insert into LR database	View online (convert to EDF/XML)	View online (via RPD)	Pass patient info to native PSG to start acquisition
Philips Respironics	Alice Sleepware G3	X	X	X	X	X	X	X	X
	Alice NightOne	X	X		X			X	
CareFusion	SomnoStar	X	X	X	X	X	X	X	X
	NOX-T3 (HST)	X	X	X	X			X	X
Natus	Grass Twin	X	X	X	X	X	X		
	Sandman	X	X	X	X	X		X	
	REMbrandt	X	X	X	X			X	
	RemLogic	X	X		X			X	X
CleveMed	Crystal PSG	X	X	X	X			X	
ResMed	Apnea Link Plus (HST)	X	X		X			X	
Cadwell	Easy III, Easy ApneaTrak	X	X	X	X			X	X
Itamar	WatchPat	X	X	X	X	X		X	X